



Philips Sonicare
HealthyWhite
Rechargeable sonic
toothbrush

3 modes

2 brush heads

HealthyWhite



HX6732/02

Whiter, healthier teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. HealthyWhite is proven to remove everyday stains in just two weeks, with the regular use of the Clean and White Mode.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Clinically proven safe and gentle

Offers you a customised cleaning experience

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2 minute mode for outstanding clean
- Sensitive Mode: For a gentle, yet thorough clean



PHILIPS
sonicare

sense and simplicity

Specifications

Cleaning performance

- Performance: Removes more plaque than Sonicare Essence/Sonicare CleanCare
- Health benefits: Helps improve gum health
- Whitening benefits: Can whiten teeth up to two shades in two weeks
- Modes: 3 Modes: Clean, Sensitive and Clean and White
- Brush heads: 2 ProResults standard brush heads for an All-Around Clean
- Speed: Up to 31,000 brush strokes per minute
- Timer: Smartimer and quadpacer

Ease of use

- Battery indicator: Illuminated icon blinks when it's time to recharge
- Brush head system: Easy snap-on brush heads for optimal hygiene
- Brushing time: 42 2-minute brushings OR 3 weeks
- Handle: Slim ergonomic design and rubber grip for easy handling.

Items included

- Additional brush head colour code rings: 3 pcs
- Brush head colour rings: 3 colour rings for easy sharing
- Charger: 1 multi-voltage charger

- Charger cover: 1 deluxe charger cover with brush head holder
- Deluxe Charger with cord wrap: 1 pcs
- Travel charger: Multi-voltage, compact travel charger

Technical specifications

- Battery: Rechargeable
- Operating time (full to empty): Full charge should provide 20 2-min. brushings

Design and finishing

- Colour: Frost White
- Translucent White front panel

Service

- Guaranteed: 2-year limited warranty

Power

- Voltage: Multi-voltage charger

Features

- ProResults brush head
- Quadpacer: Encourages attentive brushing of each quadrant of the mouth
- Smartimer: 2-minute timer helps ensure recommended brushing

Highlights

Patented sonic technology

Sonicare's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Angled brush head neck

Uniquely angled brush head neck makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Naturally whiter teeth

Dynamic Cleaning Action and direct broader contact with each tooth have been proven to remove everyday stains and help keep teeth naturally whiter.

Safe and gentle

Your Sonicare is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Smartimer

Two-minute timer automatically tracks brushing time and turns the toothbrush off after two full minutes of brushing is complete.

Quadpacer

30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Easy-start programme

Gently increases power over the first 14 uses to ease into the Sonicare experience

Clean and White Mode

2 minutes of Clean mode with an additional 30 seconds of White mode to focus on your visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks, with the regular use of the Clean and White

Clean Mode

Clean Mode: Standard 2 minute mode for outstanding clean

Sensitive Mode

Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

